

## PRODUCT SPECIFICATION SHEET

PRODUCT: Actigel ALD

PRODUCT NO.: 2701

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

<b>Analysis</b>	<b>Specification</b>
Bead Size Distribution	60 -160 $\mu\text{m}$
Flow Rate of Resin *	$\geq 6$ mL/min (200 cm/hr)
Aldehyde Groups by Titration	$\geq 50$ $\mu\text{mol/mL}$ resin

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007