

PRODUCT SPECIFICATION SHEET

PRODUCT: Actigel ALD Superflow

PRODUCT NO.: 2701S

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Bead Size Distribution	60 -160 μ m
Flow Rate of Resin *	\geq 20 mL/min (680 cm/hr)
Aldehyde Groups by Titration	\geq 50 μ mol/mL resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007