

## PRODUCT SPECIFICATION SHEET

PRODUCT: ActiClean Etox

PRODUCT NO.: 2705

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

<b>Analysis</b>	<b>Specification</b>
Endotoxin Binding Capacity	$\geq 20,000$ EU/mL resin
Bead Size Distribution	60 -160 $\mu\text{m}$
Flow Rate of Resin *	$\geq 20$ mL/min (675 cm/hr)
Amino Groups in Resin Wash	$\leq 1$ $\mu\text{mol/mL}$ resin

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 21 Aug 2008