

PRODUCT SPECIFICATION SHEET

PRODUCT: DNA Etox

PRODUCT NO.: 2706

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Endotoxin Binding Capacity	≥ 500,000 EU/mL resin
Bead Size Distribution	60 -160 μm
Flow Rate of Resin *	≥ 30 mL/min (1020 cm/hr)
BSA Binding Capacity	≤ 10 mg/mL resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007