

PRODUCT SPECIFICATION SHEET

PRODUCT: ReductEtox

PRODUCT NO.: 2709

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Endotoxin Binding Capacity	$\geq 400,000$ EU/mL resin
Bead Size Distribution	100 -200 μm
Flow Rate of Resin *	≥ 20 mL/min (675 cm/hr)
Amino Groups in Resin Wash	≤ 1 $\mu\text{mol/mL}$ resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Feb 2012