

PRODUCT SPECIFICATION SHEET

PRODUCT: Heparin Actigel

PRODUCT NO.: 37304

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Bead Size Distribution	60 -160 μ m
Flow Rate of Resin *	\geq 6 mL/min (200 cm/hr)
Factor VIII Binding Capacity	\geq 100 U/mL resin
AT3 Binding Capacity	\geq 6 mg/mL resin
Sodium Cyanoborohydride Conc.	< 0.1 μ g/mL

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 25 Nov 2013