

## PRODUCT SPECIFICATION SHEET

PRODUCT: Heparin Superflow

PRODUCT NO.: 37304S

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

<b>Analysis</b>	<b>Specification</b>
Bead Size Distribution	60 -160 $\mu\text{m}$
Flow Rate of Resin *	$\geq 20$ mL/min (680 cm/hr)
AT3 Binding Capacity	$\geq 6$ mg/mL resin
Sodium Cyanoborohydride Conc.	$< 0.1$ $\mu\text{g/mL}$

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007