

PRODUCT SPECIFICATION SHEET

PRODUCT: IDA Chelating Cellthru BigBead

PRODUCT NO.: 37311

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Bead Size Distribution	300 -500 μ m
Flow Rate of Resin *	\geq 150 mL/min (5100 cm/hr)
Hb Binding Capacity	\geq 30 mg/mL resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007