

PRODUCT SPECIFICATION SHEET

PRODUCT: Nickel Chelating Superflow

PRODUCT NO.: 37408

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

Analysis	Specification
Bead Size Distribution	60 -160 μ m
Flow Rate of Resin *	\geq 30 mL/min (1020 cm/hr)
Nickel Titration	\geq 15 μ mol/mL resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 05 Dec 2008