

PRODUCT SPECIFICATION SHEET

PRODUCT: Papain Actigel

PRODUCT NO.: 6700

STORAGE BUFFER: Acetate/Glycerol

Analysis	Specification
Bead Size Distribution	60 -160 μ m
Flow Rate of Resin *	\geq 6 mL/min (200 cm/hr)
Enzyme Activity	\geq 7 BAEE Units/mL resin
Immobilization of Papain	\geq 30 mg/mL resin

* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007