

## PRODUCT SPECIFICATION SHEET

PRODUCT: Pepsin Actigel

PRODUCT NO.: 6701

STORAGE BUFFER: Acetate/Azide/Glycerol

<b>Analysis</b>	<b>Specification</b>
Bead Size Distribution	60 -160 $\mu$ m
Flow Rate of Resin *	TBD (expect $\geq$ 6 mL/min (200 cm/hr))
Immobilization of Pepsin	$\geq$ 3 mg/mL resin

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 25 Jul 2007