

## PRODUCT SPECIFICATION SHEET

PRODUCT: Uniflow 4

PRODUCT NO.: 704

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

<b>Analysis</b>	<b>Specification</b>
Bead Size Distribution	60 -160 $\mu\text{m}$
Flow Rate of Resin *	$\geq 6$ mL/min (200 cm/hr)
Integrity of Crosslinked Beads	Pass
Size Exclusion Analysis	$^T K_{av} \geq 0.30$

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 30 May 2008