

## PRODUCT SPECIFICATION SHEET

PRODUCT: Superflow 6

PRODUCT NO.: 806

STORAGE BUFFER: 20% Ethanol (SDA 3A denatured)

<b>Analysis</b>	<b>Specification</b>
Bead Size Distribution	60 -160 $\mu$ m
Flow Rate of Resin *	$\geq$ 20 mL/min (680 cm/hr)
Macroscopic Homogeneity of Resin	Pass

---

\* Flow rate determined in a 1.5 x 5 cm column

Specification Date: 20 Jun 2007